

## ABI7300 Typical Setup Procedure

Turn on ABI7300 instrument first

Next, start the ABI7300 system software

-Window says at the bottom, “disconnected”. It will connect later.

Select File → New

-96 well

-absolute quant (std curve)

-blank document

Plate Name (e.g., Class Identifier)

Next

Select detectors (e.g., SYBR)

click Add

Passive reference = ROX

Next

Select wells that will be used and detectors in the well

Finish

System initialization should show “connected”

If not, start over

Select Plate tab

Rename wells, if you want to

-task: unknown (=experimental)

standard

NTC (no template control)

Select Instrument tab to set up reaction profile

Click on “add dissociation” (runs at very end)

Select data collection step (default should be correct already)

Save setting, or run will not occur

File → Save As... (choose your own folder, name file)

Start the run

Should take about 2 hours